



# javelin


The javelin...

is a dynamic pylon foot designed to provide excellent energy return and ground compliance. The independent heel and toe ensure a smooth and symmetrical rollover. The result is a lightweight foot for new and experienced amputees who desire a responsive foot with a sleek cosmetic profile.

Get busy living



endolite  
get busy living



Smoother than I thought possible



Let's me do more with less




Lightweight to keep up with my grandchildren



Flexible for my busy lifestyle



Keeps me in perfect balance



endolite

get busy living



# jave<sup>in</sup>

The sliding, rotating male or female adaptor, is offered separately

Dynamic pylon for tibial progression and excellent energy return

Slim and lightweight design allows for ease of cosmetic finishing

Length guide for ease of set-up. Javelin LP is cut to the minimum height prior to shipping

Independent e-carbon heel and toe springs for smooth rollover. Tripod design with split toe for maximum ground compliance



Fitting clamp allows for accurate assessment of limb performance and alignment before final cutting and drilling

Supplied fitting clamp allows the accurate fitting of the Javelin dynamic pylon to make walking relaxed and comfortable while maximizing energy return.

Initially set up the limb 10-15mm longer than required , then align to achieve perfect tibial progression, adjusting the height accordingly.

When the set-up is ideal, simply drill and cut toe blade before removing the fitting clamp. Assemble with the definitive top housing for the final fitting and limb delivery.

## Javelin Spring Selection Guide:

User weight lbs.

100	116	131	151	171	196	221	256	286	326
to	to	to	to	to	to	to	to	to	to
115	130	150	170	195	220	255	285	325	365

1 1 2 3 4 5 6 7 8 9

K3 Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

### Order example:

Javelin  
**JAV24L1M**  
size 24 left  
#1 spring kit, male  
proximal attachment  
add D suffix for dark shell



M for male attachment


Javelin Low Profile  
**JAVLP24L1F**  
size 24 left  
#1 spring kit, female proximal  
attachment  
add D suffix for dark shell



F for female attachment

### Foot Specifics:

Available Size:	22cm - 30cm
Proximal attachment	male or female
Component weight:	
*Full length	21.25 oz
*Low Profile	19.95 oz
*size 25cm without foot shell	
Max. user weight:	365 lbs
Build height:	
Full length	14.37"
Low Profile	7.87"
Heel height:	3/8"
Activity level:	3



Available with your choice of light or dark foot shell, glide sock is included. Javelin full length version ships partially assembled. Javelin LP ships fully assembled.

endolite, 1031 Byers Road, Miamisburg, OH 45342 [www.endolite.com](http://www.endolite.com)

Tel: 800.548.3534, Fax: 800.929.3636, [info@endolite.com](mailto:info@endolite.com)

Distributed by: SPS, 6025 Shiloh Road, Suite A, Alpharetta, GA 30005, Tel: 800.767.7776

[www.spsco.com](http://www.spsco.com) - <http://estore.spsco.com>