

echelon
VAC



Exceptional control and comfort



endolite

“As an amputee, socket security is one of the most important things.” Charlie

Optimal socket connection is critical to an amputee's comfort, security and stability. By pairing Biomimetic Hydraulic Technology with an elevated vacuum system, the design of EchelonVAC works to create a secure and comfortable socket connection.



Relative movement is a major issue for amputees and can lead to:

Rubbing/chafing



Potential skin breakdown

Pistoning



Increased risk of falls and reduced gait symmetry

Reduced proprioception




Less control, compromised safety

By harnessing natural ankle motion, EchelonVAC quietly creates elevated vacuum, helping to maintain an optimally fitting socket throughout the day.

With every step, the wearer presses their weight into the prosthesis, initially expelling air through a one-way valve. Simultaneously the ankle plantarflexes, actively drawing air out of the socket. This air is held in the vacuum chamber and expelled through a secondary one way valve as the tibia progresses and the ankle dorsiflexes.

The result is greater residual limb volume control and an improved connection between the residual limb and the socket. For the user, this reduces relative movement, improving proprioception and control of the prosthesis for greater comfort and safety every day.



“It’s had a great impact on my day-to-day life.” Charlie





Innovative Design

The innovative design of EchelonVAC is **lightweight** and has a **low build height** as no external power source is required.

With **no batteries** or **pump** to worry about, EchelonVAC is **quiet** and **easy to fit**.

When used in conjunction with a Silcare Breathe liner, the vacuum is applied directly to the residual limb to further enhance the connection between the limb and the socket.



The Evidence

Active vacuum systems help to stabilise residual limb volume to improve socket stability and proprioception. Scientific studies* have shown that elevated vacuum systems help to:

- Reduce volume fluctuation
- Reduce interface pressures
- Improve wound healing
- Reduce pistoning
- Improve gait symmetry, balance and reduce risk of falls
- Greater comfort and improved overall satisfaction

*Please see Clinical Compendium available to download at blatchford.co.uk/endolite/echelonvac

Features

- Biomimetic Hydraulic Technology with integrated elevated vacuum
- No power required, quiet gentle operation
- Lightweight, compact design
- Low build height
- E-carbon springs for efficient energy return
- Split toe design for ground compliance on uneven terrain
- Weatherproof - suitable for outdoor use

Technical Information

- **Maximum User Weight:** 125kg
- **Activity Level:** 2*, 3, 4*
- **Component Weight:** 930g†
- **Maximum Vacuum:** 17" Hg
- **Build Height:** Sizes 22-24: 121mm
Sizes 25-26: 126mm
Sizes 27-30: 131mm
- **Heel Height:** 10mm
- **Warranty:** 36 months
Foot shell - 12 months
Glide Sock - 3 months

Order Example

EVAC

25L 3

Size Side Spring Set


For dark tone add suffix D.
Foot example: EchelonVAC, size 25 left, spring rating 3.

Spring Set Selection

| Activity | User Weight | | | | | | | | kg lbs | Foot Spring Set |
|----------|--------------------|--------------------|--------------------|--------------------|--------------------|---------------------|----------------------|----------------------|-----------|-----------------|
| | 44-52 (100-115) | 53-59 (116-130) | 60-68 (131-150) | 69-77 (151-170) | 78-88 (171-195) | 89-100 (196-220) | 101-116 (221-255) | 117-125 (256-275) | | |
| 3 | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | | |



*Dependent on local reimbursement guidelines. May not be suitable for running or high impact activities.
†Component weight shown is for a size 26cm without footshell.

echelon
VAC



| | |
|-----------------|----------------|
| * sizes | ** sizes |
| 22 - 24 = 121mm | 22 - 24 = 70mm |
| 25 - 26 = 126mm | 25 - 30 = 75mm |
| 27 - 30 = 131mm | |

Patent numbers; US8308815, GB2536056 App, EP2124843 App, EP2124842 App, US8574312, US7985265, US8740991, US8641780, JP5336386, JP5560045, WO 2007/054736, WO 2008/071975, WO 2008/103917

| | |
|---|------------------------|
| Follow Blatchford: | blatchford.co.uk |
|  @BlatchfordGrp | 0114 263 7900 |
|  /BlatchfordGrp | sales@blatchford.co.uk |



Endolite is part of the Blatchford Group. Registered address, Chas A Blatchford & Sons Ltd, Unit D, Antura, Kingsland Business Park, Basingstoke, RG24 8PZ, United Kingdom. t: +44 (0) 1256 316600 f: +44 (0) 1256 316710



An annual visual inspection is recommended. Check for visual defects that may affect proper function. Maintenance must be carried out by competent personnel. Please check with your clinician before carrying out any new activities of daily living that you have not received training for with your EchelonVAC.