




## The Elite VT...

provides the ultimate in versatility, shock absorption, rotation and stability in one compact package.

Get busy living



endolite  
get busy living



Choletan Senior  
Fitness Instructor

"The good thing about this foot is that it allows the knee to move without the foot moving as well. I feel a lot more confidence and a lot more balance.

It's an all-rounder! Instead of having one running leg, one for rock climbing and one for walking you can have the eliteVT and do several aspects at the same time."



**eliteVT**

A precision engineered titanium spring is the core of the Elite VT's dynamic function.



The eliteVT enhances forward progression when the combined action of the axial spring and the carbon foot springs return energy at two points in the gait cycle during walking:

**Progressive transfer from heel strike to mid-stance**  
**Responsive push off**

When running this fluid action changes to a single loading response:

**Powerful loading and return of energy**



The eliteVT smoothly manages the requirements of walking and running by combining the properties of the axial spring and the foot springs. The result is a biomimetic balance similar to that provided by the combined tibial muscle and foot arch action in the natural foot.

## Elite VT Heel/Toe Selection Guide:

User weight lbs.									
100	116	131	151	171	196	221	256	286	326
to	to	to	to	to	to	to	to	to	to
115	130	150	170	195	220	255	285	325	365

**Daily Walking and Occasional Sports** such as golf and hiking

1 1 2 3 4 5 6 7 8 9

**Aggressive Walking, Frequent or Daily Sports** such as jogging

1 2 3 4 5 6 7 8 9

**Daily High Impact Activities** such as distance running, climbing, lifting, and carrying heavy objects

2 3 4 5 6 7 8 9

**Axial spring:** The Elite VT axial springs are prematched to the heel and toe spring kits by category

- Spring kits 1-3 = Axial spring 1
- Spring kits 4-5 = Axial spring 2
- Spring kits 6-7 = Axial spring 3
- Spring kits 8-9 = Axial spring 4

Standard prematched combinations will provide most active trans-tibial amputees a suitable range of vertical compression and axial rotation. You may achieve a softer or stiffer action by selecting one of the alternative springs.

### Order example:


**EVT24L21**  
size 24 left  
#2 spring kit  
#1 axial spring  
add D suffix for dark shell



Available with your choice of light or dark foot shell. Glide sock is included.  
Elite VT ships fully assembled.

### Foot Specifics:

Available Size:	24cm - 30cm
Component weight: size 26cm without foot shell	24 oz
Max. user weight:	365 lbs
Build height:	6 5/8"
Heel height:	3/8"
Activity level:	3 - 4



Endolite, 1031 Byers Rd, Miamisburg, OH 45342 [www.endolite.com](http://www.endolite.com)  
Tel: 800.548.3534, Fax: 800.929.3636, [info@endolite.com](mailto:info@endolite.com)

Distributed by: SPS, 6025 Shiloh Rd, Suite A, Alpharetta, GA 30005, Tel: 800.767.7776  
[www.spsco.com](http://www.spsco.com) - <http://estore.spsco.com>